

MARK J. STOCK
SCIENTIFIC DIGITAL ART

Framing Instructions for Show Pieces

For all 24"x24" digital C-prints

To create the standard gallery appearance, please use the following measurements and materials:

- Print (24"x24") mounted to acid-free backer board using archival mount
- 3 1/8" border, pure white 8-ply mat (8653)
- If print has a white border, include a 1/4" spacer between board and mat
- 30" x 30" frame, Nielsen 93 (matte black, square edge)
- Glazing at customer's request (plexi if piece is to be shipped)

Questions? Don't hesitate to call me (9am-11pm EST).

Mark J. Stock
734.644.8950